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OptiXcell

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Instructions for use

OptiXcell™

Please read carefully these instructions before using this product. They contain important information for the processing of bovine semen.

**OptiXcell™ is a 3 X concentrated bovine semen media used to freeze or prepare fresh bull semen.
OptiXcell™ does not contain protein of animal origin**

IDENTIFICATION OF THE PRODUCT

Composition :

Carbohydrates, mineral salts, buffer, antioxidant, phospholipids, ultra-pure water, glycerol
Antibiotics (Penicillin, Streptomycin, Lincomycin and Spectinomycin) conforming to European Regulations CEE 88/407.
OptiXcell™ does not contain protein of animal origin

Technical Characteristics

Concentrated media 3 X made according to G.M.P in the IMV Technologies media laboratory under ISO 9001 V2008.
More information on the IMV Technologies media laboratory on our website at www.imv-technologies or at the following link: <http://www.youtube.com/watch?v=EH6tjBMWJGQ>

Color: transparent beige
Product sterilized by filtration on 0.22µm filter

Presentations :

Ref : 024385 250 ml bottle. Concentrated 3 X for a final solution of 750 ml
Ref: Sample 50 ml bottle. Concentrated 3X for a final solution of 150 ml

Storage conditions

OptiXcell™ can be transported at ambient temperature (< +32°C).
For long storage, it is preferable to keep the product between +2°C and +8°C.
Do not store close to foodstuff, drinks or animal feed
Do not ingest or inhale. In case of accidental projection on skin or in the eyes, rinse profusely with water
Consult a doctor in case of problem

Expiration

18 months from the manufacturing date. (lot number and expiration date on the label).

Recommendations

The quantity of antibiotics conforms to European Regulations CEE 88/407 as long as the dilution ratio semen / extender is not over 1:10 (v:v). When doses are prepared for export countries with CSS regulations, the quantity of antibiotics has to be adjusted accordingly.

DILUTION INSTRUCTIONS:

1. Media preparation

- Place the capped bottle of OptiXcell™ in a water bath at +32/+34°C for at least 10 minutes
- Place an Erlenmeyer filled with 2 times the volume of OptiXcell™ (500 ml for a 250 ml bottle, 100 ml for a 50 ml bottle) with bi-distilled, sterile, pyrogen-free water (or equivalent quality) in a water bath at +32/+34°C for at least 10 minutes. Plug the Erlenmeyer.
- When the temperature of the water has reached +32/+34°C, remove the bottle of OptiXcell™ from the water bath. Dry it
- Slow pour the content of the OptiXcell™ bottle in the Erlenmeyer containing the ultra-pure water.
- Rinse twice the OptiXcell™ bottle off with the solution. Do not add water

The extender is now ready to use.

For a one step dilution (or if the overall process takes less than 6 hours) the entire volume can be preserved at +32/34°C.

In the case of a several steps dilution, aliquot the obtained volume into two parts. One preserved in the water bath at +32/34°C and the other part maintained at +4°C.

Note: To avoid bacterial contamination, it is imperative to use the extender for further semen dilution within 6 hours of preparation

2. Semen Pre-extension (buffering)

- After examination, delicately buffer the raw semen in the water bath (+32/+34°C) using a vol:vol (1 :1) directly into the collection tube. Incorporate the extender slowly to avoid osmotic shocks.
- Mix with a gentle swirl.
- Incubate at +32/+34°C for 10 minutes before proceeding to the final dilution.

3. Semen extension:

3.1 Single-step dilution:

- After having calculated the required extender volume, pour it in a separate container. Reserve 10 ml to rinse the collection tube off. Place the container in a water bath at +32/+34°C
- After incubation, slowly pour the semen (buffered at 1:1 as described above) into the prepared extender kept at +32/+34°C.
- Rinse the collection tube off with the 10 ml you kept on the side
Cooling methods according to type of packaging
 - ⇒ **Ambient temperature packaging**
 - ✓ Maintain the diluted semen at ambient temperature for 10 to 20 minutes before packaging
 - ✓ Package into straws at ambient temperature
 - ✓ Cool the straws off to 4°C in about 1 hour 30 minutes by placing them into a 65 mm goblet. Place a maximum of 300 mini straws (0.25) per goblet = half-filled goblet
Note : Straws placed on a freezing rack and put in a cold cabinet will reach 4°C in less than 2 minutes
 - ⇒ **Packaging at 4°C**
 - ✓ Place the flask of diluted semen inside a recipient. Add a constant level of water from the water bath + a thermometer into the recipient. Place everything into the cold cabinet at 4°C. The volume of water must be determined in order to have a progressive temperature descent to reach 4°C in about 1 hour 30 minutes. Ice can be added when the temperature has reached 12°C to accelerate the temperature descent.
Note : a cooling down using a programmable water bath can be envisioned. In this case use a slope of 1°C/min
 - ✓ Straw packaging can be done at 4°C

3.2 Multi-steps dilution

⇒ Initial dilution :

- After the initial buffering, the semen can either be diluted in a constant volume to even the descent to 4°C or in half of the required final dilution
- ✓ Place the flask of diluted semen inside a recipient. Add a constant level of water from the water bath + a thermometer into the recipient. Place everything into the cold cabinet at 4°C. The volume of water must be determined in order to have a progressive temperature descent to reach 4°C in about 1 hour 30 minutes. Ice can be added when the temperature has reached 12°C to accelerate the temperature descent.
Note : a cooling down using a programmable water bath can be envisioned. In this case use a slope of 1°C/min

⇒ Final dilution :

The balance of extender to reach the final desired concentration is added in one or several steps at 4°C. Straw packaging can then occur at 4°C.

4. Equilibration:

In all cases, the semen must stay at 4°C for 3 to 5 hours before freezing in order to stabilize the exchange between media and semen.

It is preferable to put the straws on the freezing racks just before freezing to avoid sedimentation. Or roll the straws on themselves to resuspend the cells.

5. Freezing:

Freeze the straws over LN2 vapors or in a programmable freezer such as the IMV Digitcool™, according to the following curve:

| | Temperature to reach | Slope speed (°C/mn) |
|-----------------------------------|----------------------|---------------------|
| Initial temperature of the sample | 4°C | |
| Ramp | -10°C | 5 |
| Ramp | -100°C | 40 |
| Ramp | -140°C | 20 |

For further information, please contact your local distributor or IMV Technologies at www.imv-technologies.com

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