

Duo-Scan:Go™

Rugged, smart thinking

The world's first wireless ultrasound solution designed for livestock scanning, driving improvements for better production.

Rugged, smart thinking

We are experts in making robust ultrasound scanners

IMV imaging have been pioneers in veterinary imaging for over 40 years. We design and build scanners that are rugged, reliable and built to cope with the demanding environment of the farm. Now we have developed a scanner specifically for pigs, that transmits an ultrasound image to a smart device.



Built to last

The Duo-Scan:Go is completely wireless and sealed. It has no apertures, connectors or cables. By definition, this scanner is more robust than any other scanner. This also means the Duo-Scan:Go can be easily cleaned and is bio-secure. Our technology comes with a two year warranty.

Lightweight and easy to use

Weighing just 250 grams, it is no heavier than the current probes used in the market. However, ours contains the entire ultrasound system including a rechargeable battery, a 5 hours rechargeable battery.

Reduce your scanning strain


To minimize repetitive strain injury, we have angled the 64 elements probe to make it easier to get contact with the animal.

Comfortable grip

The ergonomics of the scanner are also designed to fit everyone.

Back fat measurement

Taking back fat measurements is helpful to assess sow nutrition pre and post farrowing for better production.



Download iOS or Android app.

Contact us

Follow us

France HQ

IMV Technologies
ZI n°1 Est
61300 L'Aigle
 +33 2 33 34 64 64
 contact@imv-technologies.com
 www.imv-technologies.com

France (Angoulême)

IMV Technologies
126, Boulevard de la République,
16000 Angoulême, France
 +33 (0)5 45 92 03 57
 france@imv-technologies.com
 www.imv-technologies.com

United Kingdom

IMV Technologies
Phoenix Crescent, Strathclyde
Business Park, Scotland, ML4 3 NJ
 +44 (0)1506 460 023
 scotland@imv-technologies.com
 www.imv-technologies.com

The Netherlands

Nifa Technologies
Pallasweg 22
8938 AS Leeuwarden
 +31 58 2123279
 info@nifa.nl
 www.nifa.nl

North America

IMV Technologies USA
9501 Louisiana Ave N, Suite 300
Brooklyn Park, MN 55445
 +1 763 488 1881
 contact@imvusa.com
 www.imv-technologies.com

China

IMV Technologies China
Room 510, Building 1,
NO.800, Naxian Road,
Pudong District, Shanghai,
201203 China
 +86 21 50586918 or 50586948
 contact@imvchina.com
 www.imvchina.com

India

IMV Technologies India
(IMV India Pvt. Ltd.)
Plot No 750 Phase V, Udyog Vihar
Gurgaon 122016 (Haryana)
 +91 124 4770707
 contact@imvindia.com
 www.imvindia.com

Brazil

IMV Technologies Brazil
R. Vitoriano dos Anjos, 1081
13041-317 Campinas-SP
 +55 (19) 2513-3444
 contato@imvbrasil.com.br
 www.imv-technologies.com

IMV Technologies. For Life.
www.imv-technologies.com

