



Easi-Scan:Go[™] Lite



The Scanner You Can Rely On for Breeding Success.

Robust, versatile, and portable, this durable scanner is the little brother of our flagship product — the Easi-Scan: Go. Free from wires and built to last, the Easi-Scan:Go Lite delivers accurate images on the go — for better, more profitable breeding programs. At a lower price than its predecessors, pregnancy detection has never been more affordable.

Boost Your Breeding Programs

Perform more pregnancy detection scans per hour than ever before to boost your efficiency levels. The Easi-Scan:Go Lite is designed for speed — allowing you to cut costs and grow your profit margins. Animal pregnancy detection has never been easier — or more affordable. For high-quality ultrasound images in the most challenging environments, you can rely on the Easi-Scan:Go Lite to deliver.

Robust and Durable

Water-proof and shock-proof, this scanner is designed to withstand heavy use in busy farming environments. The detachable battery that powers the Easi-Scan:Go Lite can be removed or swapped in seconds — allowing you to scan faster than ever.

Crystal Clear Images

Capture clear, good-quality ultrasound images without lag in seconds — thanks to the Easi-Scan:Go Lite's innovative Smart Beam technology and wireless VGA BUG goggles.

**UP TO 3 SMART DEVICES
OR WIRELESS VGA BUG**
can be connected

WATERPROOF
Sealed buttons with no
connectors

850 GRAMS

IMPROVED FLEXIBILITY
for better cable
management



USE CALIPERS
to measure the embryo
size

RECORD VIDEOS
and save images
on your smart devices

**LOW COST LITHIUM-ION
BATTERY**
up to 5 hours

128 ELEMENTS
receiving on 8
channels, 26 frames
per seconds. Linear
and Curved probes
available.

Easi-Scan:Go Lite Unique Technology



Enhanced Efficiency

Sound-Link IMV's unique wireless data transmission platform for real time imaging allows image display without any image lag. With Smart-Beam technology we achieve the lowest power ultrasound scanner with standard Wi-Fi devices used as a primary viewing device. Sharing images on multiple devices allows for better engagement with farm staff.



Unrivaled Longevity

The tough farm conditions require reliable equipment, which is why we developed Shield-Tech, that makes the Easi-Scan:Go Lite our most rugged scanner yet. Our fixed probe with strengthened cable has improved durability. The scanner can be washed in a bucket of water improving workflow and biosecurity.



Productivity

Used in conjunction with the new Visits ToDo App, the Easi-Scan:Go Lite captures cow-side data, making clipboard recording a thing of the past. With visual prompts for both the vet and farmer you can ensure diagnostic results are recorded correctly.

Contact us

Follow us   

France HQ

 IMV Technologies
ZI n°1 Est
61300 L'Aigle
 +33 2 33 34 64 64
 contact@imv-technologies.com
 www.imv-technologies.com

France (Angoulême)

 IMV Technologies
126, Boulevard de la République,
16000 Angoulême, France
 +33 (0)5 45 92 03 57
 france@imv-technologies.com
 www.imv-technologies.com

United Kingdom

 IMV Technologies
Phoenix Crescent, Strathclyde
Business Park, Scotland, ML4 3 NJ
 +44 (0)1506 460 023
 scotland@imv-technologies.com
 www.imv-technologies.com

The Netherlands

 Nifa Technologies
Pallasweg 22
8938 AS Leeuwarden
 +31 58 2123279
 info@nifa.nl
 www.nifa.nl

North America

 IMV Technologies USA
9501 Louisiana Ave N, Suite 300
Brooklyn Park, MN 55445
 +1 763 488 1881
 contact@imvusa.com
 www.imv-technologies.com

China

 IMV Technologies China
Room 510, Building 1,
NO.800, Naxian Road,
Pudong District, Shanghai,
201203 China
 +86 21 50586918 or 50586948
 contact@imvchina.com
 www.imvchina.com

India

 IMV Technologies India
(IMV India Pvt. Ltd.)
Plot No 750 Phase V, Udyog Vihar
Gurgaon 122016 (Haryana)
 +91 124 4770707
 contact@imvindia.com
 www.imvindia.com

Brazil

 IMV Technologies Brazil
R. Vitoriano dos Anjos, 1081
13041-317 Campinas-SP
 +55 (19) 2513-3444
 contato@imvbrasil.com.br
 www.imv-technologies.com