Gallicomfort

Turkey insemination bench (all options included)	ref. 027988
Turkey insemination bench (without options)	ref. 027983

Options

1 MRA holder	ref. 027984
1 ergonomic arm rest for the inseminator	ref. 027985
1 release ramp	ref. 027986
1 stool for the inseminator	ref. 027987

Related products

2.5ml poultry media - antibiotic free	ref. 002262
4ml poultry media - antibiotic free	ref. 020070
5ml Extendyl	ref. 025880
70mm straws	ref. 006929
MRA filling machine	ref. 007155
Motorized kit	ref. 019250
Electric block	ref. 008593
Pump syringe for MRA filling machine	ref. 007161
Poultry semen processing workshop	ref. 024105

Gallicomfort

Turkey insemination bench









What is Gallicomfort?

Gallicomfort is a turkey insemination bench designed by an ergonomic specialist to help improve both human and animal welfare, while maintaining productivity levels.





- 1 Turkeys are guided up an access ramp to the insemination bench.
- **2** Breeding crew members transfer hens, one at a time, at lifting height, to the inseminator by sliding them across the bench's table.
- **3** The inseminator, sitting on a stool, controls the hen by pressing the foot pedal. This engages the leg holder and rollers to open the hen's oviduct.
- **4** The ergonomic arm rest allows the inseminator to easily guide his arm from the MRA to the hen with a filled semen straw for insemination.
- **5** The inseminator then presses the pedal again to release the hen. A ramp positioned below the bench allows for a gentle release.

Why choose Gallicomfort?

Put your staff's well-being first

Gallicomfort is designed to improve your staff's working conditions and reduce the risk of musculoskeletal disorders:

- Turkeys arrive at lifting height (no bending)
- Carry weight reduced by half
- Ergonomic arm rest
- Adjustable stool



Classical AI AI with Gallicomfort

Study carried out on a turkey farm in Brittany, for 2,400 inseminations



Keep your efficiency in mind

Maintain high productivity levels

• 700 turkeys / hour with a team of 5 people

Train staff quickly

• Train new employees in 1 to 2 weeks

Make animal welfare a priority

Gallicomfort is designed to improve both employees and animal welfare.

- Reduce the use of small pens, allowing for hens to move around freely.
- Gentle handling: turkeys are softly slid onto the bench for insemination. They are released calmly down a ramp after insemination.

