

Connected to  
... PROTECTION

# Shoulder Protect glove

FOR FULL PROTECTION




- Full shoulder, neck and armpit protection
- Ideal for exploration and search
- Highly resistant




# Shoulder Protect glove

FOR FULL PROTECTION



## The complete shoulder and armpit protection

- Neck strap holds glove firmly in place
- Protects against excessive soiling without gown use
- Ideal for veterinary acts (distocia, abortion, C section, prolapsus, metritis, ultrasound...) or multiple AI's and palpations



### Ref. 025593 – Shoulder protect glove

Polyethylene gloves

Thickness: 30 microns 1.2 mil ''

Clear color

Man size

Box of 50 gloves

13 x 50 gloves per box