### **Connected to**

... AQUACULTURE REPRODUCTION



# Aquaculture CryoLab

## READY-TO-USE SPERM ALL-YEAR ROUND

- Fast, reliable machines and procedures to freeze fish milt
- Fertilization of 300 eggs with one straw (0.5ml)
- Semen identification and traceability
- Up to 50 milion eggs fertilized with milt stored in one tank
- 90 % fertilization rate with straws







## Aquaculture CryoLab

READY-TO-USE SPERM ALL-YEAR ROUND

#### Packing and sealing

- Automated detection of plug saturation
- Optimum filling from the first to the last straw
- Minimized human intervention



ISevo



#### Automated freezers

- • • • • • • • • •
- State of the art solution based on sealed cryogenic containers (-150°C to +50°C)
- Automated recording of freezing curves for each cycle
- Repeatable freezing cycles
- Easy freezing curve parameter definition
- Four sizes: 33, 240, 1392, 3000 straws

#### Digitcool

#### Long-term storage

- 400 liters container
- 145 600 straws in one tank
- Identification of goblets and straws
- Capacity to fertilize 2 000 4 000 females



CF 400 nitrogen tank

www.imv-technologies.com