

Duo-Scan:Go™



Rugged, smart thinking

The world's first wireless ultrasound solution designed for livestock scanning, driving improvements for better production.

Rugged, smart thinking

We are experts in making robust ultrasound scanners

IMV imaging have been pioneers in veterinary imaging for over 40 years. We design and build scanners that are rugged, reliable and built to cope with the demanding environment of the farm. Now we have developed a scanner specifically for pigs, that transmits an ultrasound image to a smart device.



Built to last

The Duo-Scan:Go is completely wireless and sealed.
It has no apertures, connectors or cables. By definition, this scanner is more robust than any other scanner. This also means the Duo-Scan:Go can be easily cleaned and is bio-secure.
Our technology comes with a two year warranty.

Lightweight and easy to use

Weighing just 250 grams, it is no heavier than the current probes used in the market. However, ours contains the entire ultrasound system including 5 hours rechargeable battery.

Reduce your scanning strain


To minimize repetitive strain injury, we have angled the 64 elements probe to make it easier to get contact with the animal.

Comfortable grip

The ergonomics of the scanner are also designed to fit everyone.

Back fat measurement

Taking back fat measurements is helpful to assess sow nutrition pre and post farrowing for better production.



Go

Download iOS or Android app.

Contact us

Follow us

France

IMV Technologies Headquarters
+33 2 33 34 64 64
61300 - L'Aigle
welcome@imv-technologies.com

IMV imaging
16000 - Angoulême
france@imv-technologies.com

Polysem
78120 - Rambouillet
welcome@imv-technologies.com

The Netherlands

Leja
Nieuw-Vennep - 2153 GN - North Holland
info@leja.nl

NIFA Technologies
Leeuwarden - 8938 AS - Friesland
www.nifa.nl

Scotland

IMV Technologies UK
Bellshill - ML4 3NJ - North Lanarkshire

Spain

Imporvet
Barcelona - 08540 Centelles - Catalonia
int@import-vet.com

United States

IMV Technologies USA
Brooklyn Park - 55445 - Minnesota
contact@imvusa.com

IMV Technologies USA
Rochester - 55901 - Minnesota
info@imvusa.com

Canada

Ro-main
Saint-Lambert-de-Lauzon
G0S 2W0 - Quebec
info@ro-main.com

Brazil

IMV Technologies Brazil
Campinas - 13030-300 - Sao Paulo
contato@imvbrasil.com.br

China

IMV Technologies China
201203 - Shanghai
contact@imvchina.com

Nelson Techno Medical Co
Changzhou - Jiangsu - China
contact@imvchina.com

India

IMV Technologies India
Gurugram - 122016, Haryana
contact@imvindia.com

www.imv-technologies.com

Ref: