

IMAGo™

Designed to keep you on the go

imv TECHNOLOGIES

IMAGo™

Introducing the new generation of ImaGo fertility and pregnancy scanners.

Offering unparallel accuracy and ease of use.


Designed to keep you on the go

ImaGo is tailored for bovine vets on the go. ImaGo's rugged, portable design and good image quality ensures you're always ready for high-quality care, when working across multiple farms and to busy schedules thus reducing the time and cost up to 10%.



Quick start and seamless scanning process

- Screen form factor and a choice of the pre-sets species-specific configurations ensures quick start process that only takes couple of minutes to start.



- Comfortable: wrist or chest mount options in combination with lightweight (1.18 kg) allows to reduce strain during long days and speed the scanning process. No need to worry about recording the time on the farm, ImaGo does it for you. It switches on when you start scanning automatically.





- Versatile probes: 128 elements curved or linear probes versions to suit different needs.



imv TECHNOLOGIES

- Sound-Link technology ensures connection with mobile device and goggles to share or save the ultrasound image, train students or engage with your customers.




Built to keep you scanning



- ImaGo's long-lasting, interchangeable battery provides up to four hours of continuous use and can be charged on-the-go in your car between farms due to the Smart-Beam technology.

- Shield-Tech technology makes ImaGo fully waterproof, designed to IP67 standard and rugged - eliminating downtime caused by damage, or wear-and-tear. It ensures biosecurity and saving the time of cleaning when moving between farms.





Want to detect pregnancies earlier?

- With colour flow doppler mode ImaGo helps detect non-pregnant cows earlier savin up to £3 per day per cow.

Efficiency you can rely on

- ImaGo's simple interface and innovative features help you **get more done with less effort**. Speed up your workflow, reduce downtime, and stay flexible, no matter where your day takes you.

Testimonial

You can immerse ImaGo scanner fully in water to clean it, which really helps with our bio security and making sure that we're not transferring any unwanted pathogens from farm to farm.

Victoria White,
Veterinary Surgeons
FarmVets South West, UK



imv TECHNOLOGIES

IMV TECHNOLOGIES. FOR LIFE.



- Our IMV Technologies offices
- Main distributors

Our global network: more than 120 distributors supporting you worldwide

France

IMV Technologies Headquarters
+33 2 33 34 64 64
61300 - L'Aigle
welcome@imv-technologies.com

IMV imaging
16000 - Angoulême
france@imv-technologies.com

Polysem
78120 - Rambouillet
welcome@imv-technologies.com

The Netherlands

Leja
Nieuw-Vennep - 2153 GN - North Holland
info@leja.nl

NIFA Technologies
Leeuwarden - 8938 AS - Friesland
www.nifa.nl

Scotland

IMV Technologies UK
Bellshill - ML4 3NJ - North Lanarkshire

Spain

Imporvet
Barcelona - 08540 Centelles - Catalonia
int@import-vet.com

United States

IMV Technologies USA
Brooklyn Park - 55445 - Minnesota
contact@imvusa.com

IMV Technologies USA
Rochester - 55901 - Minnesota
info@imvusa.com

Canada

Ro-main
Saint-Lambert-de-Lauzon
G0S 2W0 - Quebec
info@ro-main.com

Brazil

IMV Technologies Brazil
Campinas - 13030-300 - Sao Paulo
contato@imvbrasil.com.br

China

IMV Technologies China
201203 - Shanghai
contact@imvchina.com

Nelson Techno Medical Co
Changzhou - Jiangsu - China
contact@imvchina.com

India

IMV Technologies India
Gurugram - 122016, Haryana
contact@imvindia.com

www.imv-technologies.com

Ref: